



- “Off-the-shelf” availability for craniofacial reconstruction and augmentation
- Saves the time and expense associated with harvesting graft material
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 year shelf life



CAT#	A	B	C
4001	38mm	50mm	0.25mm
4002	50mm	76mm	0.25mm
4003	38mm	50mm	0.35mm
4004	50mm	76mm	0.35mm
4005	30mm	50mm	0.40mm
4006	38mm	50mm	0.45mm
4007	50mm	76mm	0.45mm
4008	38mm	50mm	0.85mm
4009	50mm	76mm	0.85mm
4012	38mm	50mm	1.5mm
4013	50mm	76mm	1.5mm
4015	38mm	50mm	3.0mm

