


- Excellent option for craniofacial reconstruction and augmentation
- Available in multiple sizes and thicknesses to ensure the surgeon has a proper amount of biomaterial for almost any possible contour or shape required to fit the need of the individual patient
- May be carved using a scalpel, burr, or cut with scissors
- Screws or sutures can be placed in any desired location
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	A	B	C
4019	13mm	38mm	3mm
4020	29mm	50mm	3mm
4021	38mm	63mm	3mm
4022	13mm	38mm	6mm
4023	25mm	50mm	6mm
4034	38mm	36mm	6mm
4025	13mm	38mm	9.5mm
4026	25mm	50mm	9.5mm
4027	38mm	63mm	9.5mm

