

- Restores volume of the orbit following enucleation and evisceration procedures
- Pre-fabricated tunnels for suture to pass through easily
- Tunnels allow for easy rectus muscle attachment without the use of an implant wrap
- Features a seamless transition from open porous body to a smooth anterior cap to minimize late implant exposures
- Improved tunnel design with internal ridges to assist in needle support
- Available in Spheres and Conical Orbital Implants (COI)
- COI provides approximately equivalent volume of a sphere with 2mm larger diameter
- Enucleation Introducer provided sterile with each implant
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



Sphere

CAT#	Diameter
4036	16mm
4037	18mm
4039	20mm
4041	22mm

Conical Orbital Implant (COI)

CAT#	Volume	Diameter
4059	3.0 mL	16mm
4060	4.2 mL	18mm
4062	5.6 mL	20mm
4063	7.4 mL	22mm

