

- Restores volume of the orbit following enucleation and evisceration procedures
- Available in sphere and conical shapes
- Conical Orbital Implant (COI) provides more volume posteriorly
- COI provides approximately equivalent volume of a sphere with 2mm larger diameter
- Su-Por biomaterial offers superior, unmatched suturability
- Fully porous sphere promotes tissue ingrowth
- Wrapping of implant is optional, as extraocular muscles may be sutured directly to sphere
- Enucleation Introducer provided sterile with each implant
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



Conical Orbital Implant (COI)		
CAT#	Volume	Diameter
4054	3.0 mL	16mm
4055	4.2 mL	18mm
4057	5.6 mL	20mm
4058	7.4 mL	22mm

Sphere	
CAT#	Diameter
4028	14mm
4029	16mm
4030	18mm
4031	19mm
4032	20mm
4033	21mm
4034	22mm
4035	23mm

