


- Provides up to 5mm of anterior projection
- Small flange mirrors the contour of the most anterior region of the orbital floor allowing for proper positioning
- Designed to be trimmed to meet the needs of the patient
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4064	Left	43mm	18mm	3.2mm
4065	Right	43mm	18mm	3.2mm

