


## Extended Orbital Rim

Designed with Robert A. Goldberg, M.D.

- Fits a large array of orbital rim defects
- Provides surgeons with a large area of implant
- Designed to be trimmed to meet the needs of the individual patient
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4066	Left	47mm	40mm	6.33mm
4067	Right	47mm	40mm	6.33mm

