


SP1 Malar

Designed with Oscar M. Ramirez, M.D.

- Subtly augments the malar bone while accommodating the inferior orbital nerve
- Does not significantly increase the bizygomatic distance
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4082	Small - Left	50mm	19mm	3mm
4083	Small - Right	50mm	19mm	3mm
4084	Regular - Left	50mm	19mm	5mm
4085	Regular - Right	50mm	19mm	5mm

