



- Subtly augments the malar bone
- Extends from the zygomatic arch to the maxillary buttress, below the inferior orbital nerve
- Can be used for augmentation of the malar or for traumatic reconstruction of the zygoma
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4086	Small - Left	64mm	19mm	3mm
4087	Small - Right	64mm	19mm	3mm
4197	Medium - Left	64mm	19mm	4.5mm
4198	Medium - Right	64mm	19mm	4.5mm
4199	Large - Left	64mm	19mm	7mm
4200	Large - Right	64mm	19mm	7mm

