


- For small to medium sized cranial defects and deformities
- Designed with a series of pedicles on the interior surface to provide volume and flexibility
- Available with or without non-integrating membrane layer
- Screws or sutures can be placed in any desired location
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4109	Regular	56mm	91mm	4.5mm
4134	Regular + Membrane	56mm	91mm	4.5mm

