


- Augments deficient soft tissue in the temporal region
- Thin contoured temporal surface of Temporal Flex Grid is designed for proper anatomical fit and more natural result
- SP Temporal Flex Grid designed with flat temporal surface for universal fit
- Screws or sutures can be placed in any desired location
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



SP Temporal Flex Grid



Temporal Flex Grid



| CAT# | Description    | A    | B     | C    |
|------|----------------|------|-------|------|
| 4110 | Small - Left   | 61mm | 78mm  | 18mm |
| 4111 | Small - Right  | 61mm | 78mm  | 18mm |
| 4112 | Medium - Left  | 74mm | 93mm  | 20mm |
| 4113 | Medium - Right | 74mm | 93mm  | 20mm |
| 4114 | Large - Left   | 82mm | 105mm | 20mm |
| 4115 | Large - Right  | 82mm | 105mm | 20mm |
| 4363 | SP- Small      | 70mm | 70mm  | 10mm |
| 4364 | SP -Medium     | 88mm | 86mm  | 15mm |
| 4365 | SP - Large     | 98mm | 95mm  | 18mm |

