


Burr Hole Cover


- Designed to fit into and over holes made by a cranial perforator
- Large Burr Hole Cover stem is designed for 14mm diameter holes and can easily be modified to fit smaller holes
- Small Burr Hole Cover stem designed for 5mm diameter holes
- Screws or sutures can be placed in any desired location
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- 3 per package, MR Safe, Sterile, 5 Year Shelf Life



Small



Large



CAT#	Description	A	B	C
4130	Small	15mm	3mm	5mm
4131	Large	29mm	7mm	14mm

