



## Enophthalmos Wedge

- Mimics the shape and contour of the orbital floor
- Restores the volume and shape of the orbit
- Fully porous to promote tissue ingrowth and facilitate implant stabilization
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4180	Regular - Left	24mm	33.5mm	7mm
4181	Regular - Right	24mm	33.5mm	7mm
4182	Large - Left	28mm	40mm	7.5mm
4183	Large - Right	28mm	40mm	7.5mm

