


Midface Rim

- Augments or repairs the non-load bearing bony structures of the midface and inferior orbital rim
- May be used in its entirety or trimmed to meet the needs of the individual patient
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4189	Left	47m	28mm	3mm
4190	Right	47m	28mm	3mm

