


- “Off-the-shelf” availability for orbital floor reconstruction and augmentation
- Save the time and expense associated with harvesting graft material
- Composite structure comprised of a porous bone-side layer heat-bonded to a non-porous, non-integrating plate
- Fully porous bone-side layer promotes tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	A	B	C
4231	38mm	50mm	0.6mm
4016	38mm	50mm	1.0mm
4241	38mm	50mm	1.6mm
4017	50mm	76mm	1.0mm
4018	50mm	76mm	1.6mm

