



Flor-Tec[™] Composite Orbital Floor

Designed with William R. Nunery, M.D., F.A.C.S.
U.S. Patent No. US09724198

- “Off-the-shelf” availability for orbital floor reconstruction and augmentation
- Save the time and expense associated with harvesting graft material
- Composite structure comprised of non-porous primary section and leading porous strip
- Fully porous strip promotes tissue ingrowth and facilitates implant stabilization
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	A	B	C
4233	38mm	50mm	0.3mm
4232	38mm	50mm	0.4mm
4234	38mm	50mm	0.5mm

