

Nasal Dorsum

Petite Nasal Dorsum Designed with:
Randal Tanh Hoang Pham, M.D, F.S., F.A.C.S.


- Subtle augments the dorsum
- Tapered profile provides added height to dorsum
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life




Petite



Small / Large



Small / Large



Petite

CAT#	Description	A	B	C
4239	Small	54mm	6mm	11mm
4238	Large	67mm	9mm	14mm
4313	Petite	4mm	4mm	45mm
4314	Petite	4mm	4mm	55mm
4315	Petite	5mm	5mm	45mm
4316	Petite	5mm	5mm	55mm
4317	Petite	6mm	9mm	55mm

