


Nasal Arch

- Designed to create a nasal onlay for subtle dorsum augmentation
- Must be placed appropriately in the dorsum area and avoid extending proximally into the soft nasal cartilage area of the tip
- Edges should be feathered for a smooth transition
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually Packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4244	Small	70mm	13mm	2mm
4245	Medium	70mm	15mm	2mm
4246	Large	70mm	17mm	2mm

