


Orbital Rim Onlay

- Augments lateral and inferior orbital rims
- Subtly increases anterior projection
- Can be used in its entirety or trimmed to meet the needs of the individual patient
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4253	Left	40mm	40mm	8.45mm
4254	Right	40mm	40mm	8.45mm

