



- For significant orbital floor and wall trauma
- Addition of rigid fixation plates provide increased structural support
- Available with or without heat-bonded non-porous membrane layer
- Flat posterior side can be rested along the posterior border of the defect
- Fully porous bone-side layer promotes tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



| CAT# | Description                       | A    | B    | C     | Channel Width |
|------|-----------------------------------|------|------|-------|---------------|
| 4298 | Miniplate Channel Sheet           | 40mm | 52mm | 2.3mm | 4.0mm         |
| 4299 | Microplate Channel Sheet          | 40mm | 52mm | 2.3mm | 2.3mm         |
| 4300 | Membrane Miniplate Channel Sheet  | 40mm | 52mm | 2.3mm | 4.0mm         |
| 4301 | Membrane Microplate Channel Sheet | 40mm | 52mm | 2.3mm | 2.3mm         |

