




- Excellent reconstruction option for correcting nasal deformities
- Mimics the shape of the nasal bones and upper lateral cartilage
- Nasal Shell inserts can be used in dorsal areas where additional augmentation is required
- Identical sterile implant included for use as a template
- Sutures may be placed in any desired location
- Easy to trim with scissors or scalpel
- MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4368	Regular	38mm	21mm	17mm
	Insert - Small (included)	30mm	9mm	4mm
	Insert - Large (included)	38mm	9mm	2.5mm
4369	Large	40mm	20mm	18mm
	Insert - Small (included)	32mm	9mm	4mm
	Insert - Large (included)	41mm	9mm	3.1mm

