


Lower Eyelid Spacer

Designed with Michael A. Burnstine, M.D.

- “Excellent addition to the management of patients with symptomatic paralytic ectropion” - *Michael A. Burnstine, M.D.*
- Thin, tissue-like profile for greater eyelid stability
- Fully porous to promote tissue ingrowth and facilitate implant stabilization
- Placed inferior to the tarsus under no tension
- Available with strip to provide additional support
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4436	Lower Eyelid Spacer	29.5mm	9mm	0.45mm
4437	Lower Eyelid Spacer w Superior Strip	29.5mm	9mm	0.45mm

Scan for Video
Demonstration

