


Capocelli Plate

Designed with Anthony L. Capocelli, M.D.

- Restores cranium following cranial decompression procedures
- Screws or sutures can be placed in any desired location
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Two per package, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4461	Small	43mm	36mm	3mm
4462	Medium	48mm	40mm	3mm
4463	Large	54mm	45mm	3mm

