

Patient-Specific Chin Implant

- Augment or Restore Volume for a Natural Aesthetic Appearance
- Designed using Patient CT Data for a Precise Fit
- Designed to Meet Surgeon and Patient Needs (preoperative planning provided for all cases, including bone, tumor, and/or hardware removal)
- Fully Porous Material Facilitating Tissue Ingrowth
- Easily Modified with Scissors, Scalpel, or Burr In Vivo
- Individually Packaged, MR Safe, Sterile
- 5 Year Shelf Life
- Mandibular Implant not Included

*Leader in Quality.
Leader in Value.
Individually Crafted.*

