


SU-POR[®]


Surgical Implants

Pterional

- Corrects temporal hollowing in patients following a Pterional approach to the brain
- Place deep in the temporalis during closure
- Fully porous material for tissue ingrowth
- Available in Left /Right configurations, Smooth – Left/Right Configurations
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf life



Right Side Shown
SMOOTH



CAT #	Description	A	B	C
4122	Left	6mm	43mm	44mm
4123	Right	6mm	43mm	44mm
4467	Smooth - Left	46.3mm	43mm	2.25mm
4468	Smooth - Right	46.3mm	43mm	2.25mm



Toll Free: (877) 631-1954
www.su-por-ior.com
 © Poriferous, LLC PSS4122 REV B

PORIFEROUS[®]