




# Mastoid

- Provides an excellent option for repair of mastoid defects
- Regular Mastoid available in left and right configurations
- Small / Large Mastoid provides universal left / right fit
- Designed to be trimmed to meet the needs of the individual patient and overlap the defect
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Screws or sutures can be placed in any desired location
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4132	Small	36mm	45mm	1.5mm
4133	Large	55.2mm	63mm	1.0mm
4124	Left	56mm	53.3mm	1.5mm
4125	Right	56mm	53.3mm	1.5mm

