



Single Channel Sheets

- For significant orbital floor and wall trauma
- Addition of single rigid fixation plate provides increased structural support
- Available with or without heat bonded non-porous membrane layer
- Flat posterior side can be rested along the posterior border of the defect
- Fully porous bone-side layer promotes tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C	Channel Width
4297	Microplate Single Channel Sheet	38mm	50mm	0.85mm	2.3mm
4296	Membrane Microplate Single Channel Sheet	38mm	50mm	0.85mm	2.3mm

