



Retrosigmoid

Designed with Ali R. Zomorodi, M.D.

- The Retrosigmoid is an excellent option for use after a suboccipital craniotomy. It is available in left/right, regular/large, and low profile configuration combinations.
- The Retrosigmoid Regular is designed for typical approach.
- The Retrosigmoid Large is designed for the extended approach.
- Both configurations incorporate an impervious layer bonded to the dura facing side of the implant.
- Typical fixation utilizes 1.5mm self-drilling screws, (not included) saving time in the OR
- Retained Shape Memory (RSM)[™] provides malleability by hand.



CAT #	Description	A	B	C
4562	Low Profile - Right	29mm	35mm	.85mm
4563	Low Profile - Left	29mm	35mm	.85mm
4564	Low Profile - Large - Right	52.6mm	58mm	.85mm
4565	Low Profile - Large - Left	52.6mm	58mm	.85mm
4566	Right	32.2mm	35mm	1.6mm
4567	Left	32.2mm	35mm	1.6mm
4568	Large - Right	52.6mm	58mm	1.6mm
4569	Large - Left	52.6mm	58mm	1.6mm

