



3D Pterional

- The Pterional Implant is designed to correct temporal hollowing
- The implant is placed deep to the temporalis during closure to correct this defect.
- Sutures easily pass through the implant allowing for tissue suspension.
- Fixation holes for use with 1.2mm x 5mm long self-drilling screws.
- Shape Retention Memory[™] provides malleability by hand.
- Provided sterile by ETO
- MR Safe, 5 Year Shelf Life



Right Side Shown

Cat	Description	A	B	C
4469	3D Pterional - Left	44mm	48mm	6mm
4470	3D Pterional - Right	44mm	48mm	6mm

