


Ear Wedge

Designed with Arturo Bonilla, M.D.

- Available in both full-thickness and thin structure
- Thin structure provides improved fit to the postauricular sulcus
- Enhances the projection of a surgically reconstructed ear
- Can be trimmed to match the projection of the contralateral ear
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT #	Description	A	B	C
4292	Right	15mm	44mm	11mm
4293	Left	15mm	44mm	11mm
4294	Thin - Left	15mm	44mm	11mm
4295	Thin-Right	15mm	44mm	11mm
4285	Extended Thin - Regular - Left	15mm	44mm	11mm
4286	Extended Thin - Regular - Right	15mm	44mm	11mm
4287	Extended Thin - Large - Left	15mm	53mm	11mm
4288	Extended Thin - Large - Right	15mm	53mm	11mm

