



### Surgical Implants

- Augments the ramus of the mandible for significant augmentation
- Unmatched shape retention and suturability
- Fully porous material for tissue ingrowth
- Screws or sutures can be placed in any desired location
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C	D	E
4095	Large 5mm - Left	79mm	32mm	5mm	10mm	7mm
4096	Large 5mm - Right	79mm	32mm	5mm	10mm	7mm
4097	Large 10mm - Left	79mm	32mm	10mm	10mm	7mm
4098	Large 10mm - Right	79mm	32mm	10mm	10mm	7mm
4145	Regular 5mm - Left	79mm	32mm	5mm	5mm	4mm
4145	Regular 5mm - Right	79mm	32mm	5mm	5mm	4mm
4147	Regular 10mm - Left	79mm	32mm	10mm	5mm	4mm
4148	Regular 10mm - Right	79mm	32mm	10mm	5mm	4mm

