




Quadro-Port Tunnel Sphere

Surgical Implants

Designed with Steven C. Dresner, M.D.

- Restores volume of the orbit following enucleation and evisceration procedures
- Pre-fabricated tunnels for suture to pass through for easy attachment of extraocular muscles
- Smooth anterior surface
- Introducer provided sterile with each implant
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Diameter
4036	16mm
4037	18mm
4039	20mm
4041	22mm

Toll Free: (877) 631-1954
www.su-por-ior.com

© Poriferous, LLC PSS4036REV B

PORIFEROUS®