


Surgical Implants

USA & International Patent Protected

- Restores volume of the orbit following enucleation and evisceration procedures
- A smooth, suturable, porous exterior encasing a non-integrating orbital sphere
- Tissue ingrowth potential is limited to the thin porous outer layer eliminates the potential for a center space void of fibrovascular ingrowth
- Wrapping of implant is optional, as extraocular muscles may be sutured directly to sphere
- Introducer provided sterile with each implant
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



Introducer



Shown in half for clarity



CAT#	Diameter
4042	16mm
4043	18mm
4045	20mm
4047	22mm