



### Surgical Implants

- Subtly augments the malar bone
- Extends from the zygomatic arch to the maxillary buttress, below the inferior orbital nerve
- Can be used for augmentation of the malar or for traumatic reconstruction of the zygoma
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Screws or sutures can be placed in any desired location
- Easy to trim with fluted burr or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



| CAT# | Description    | A    | B    | C     |
|------|----------------|------|------|-------|
| 4086 | Small - Left   | 64mm | 19mm | 3mm   |
| 4087 | Small - Right  | 64mm | 19mm | 3mm   |
| 4197 | Medium - Left  | 64mm | 19mm | 4.5mm |
| 4198 | Medium - Right | 64mm | 19mm | 4.5mm |
| 4199 | Large - Left   | 64mm | 19mm | 7mm   |
| 4200 | Large - Right  | 64mm | 19mm | 7mm   |

