


Surgical Implants

- For significant orbital floor and wall trauma
- Addition of rigid fixation plates provide increased structural support
- Available with or without heat-bonded non-porous membrane layer
- Flat posterior side can be rested along the posterior border of the defect
- Fully porous bone-side layer promotes tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



Plates Shown. Not Included



Channel Details



Micro and Mini-plate sized channels available

CAT#	Description	A	C
4298	Miniplate - Regular	40mm	52mm
4299	Microplate - Regular	40mm	52mm
4300	Miniplate - Membrane	40mm	52mm
4301	Microplate - Membrane	40mm	52mm

